

Extending a Scintillation Counter's Dynamic Range

Lewis Carroll

Carroll & Ramsey Associates
Berkeley, CA, USA

Introduction Our compact, solid-state scintillation probes are widely used as HPLC / GC radiation detectors for quality assurance in PET/nuclear medicine research labs and radio-pharmacies. The detector probes operate in AC-coupled, *pulse-counting mode*, with a threshold discriminator to exclude noise and to minimize baseline fluctuation and drift.

The threshold discriminator is followed by an analog ratemeter to produce a voltage signal that is proportional to the time-rate of photon-induced pulses which exceed the pre-set threshold. Using this scheme, the ability to discern and evaluate the smallest radio-chromatography peaks – the minimum detectable signal – is governed by fluctuations in the base-line from ambient radiation background in the lab which, in turn, requires that the detector probe be well shielded so that it 'sees' only the radiation emanating from a loop of flow-tubing placed in tight proximity to the probe.

While this scheme is optimum for detection at low-to-moderate levels of radioactivity encountered in a typical quality-assurance radio-assay, pulse-counting detectors generally suffer from saturation effects due to counting system *dead-time* when exposed to high levels of radioactivity. In an effort to broaden the potential application of our scintillation detector products, we are engaged in an ongoing development program to enhance detector system linearity and dynamic range by reducing saturation effects at the 'high-end' while preserving system sensitivity at the 'low end'.

Stress-Testing at high count-rates To facilitate our development, we use home-made random pulse generators ¹ operating in parallel. Each pulse generator drives its own light-emitting diode to simulate scintillation pulses (pulse width ~ 200 nsec) from a CsI(Tl) scintillator crystal. The fixed-amplitude, random light-pulses are pre-set to match the 511 KeV principal peak in our 1 cm³ crystal, and are directed at a 1 cm² Si PIN diode + charge-integrating preamplifier (to include the effects of electronic noise inherent in a room-temperature semiconductor diode detector) all placed inside a light-tight enclosure to emulate our scintillation detector probe's 'front end'. Each generator delivers pulses at Poisson random intervals with an adjustable mean rate covering a range of ~100 pulses per second up to ~125K pulses per second. A pair of generators can produce a mean rate up to ~250K pulses per second, providing a convenient, readily-controllable source of detector system excitation over a wide range of count-rates, without having to handle large quantities of radioactive material. The 'Poisson-ness' of our random pulse generators was validated by recording the distribution of inter-pulse waiting times for various mean rates, using a calibrated time-to-amplitude converter plus multi-channel analyzer.

Extending Dynamic range In a radiation *counter*, input pulses which exceed a pre-determined threshold generate corresponding output pulses of fixed amplitude which, in turn, are either counted digitally or time-averaged in an analog rate-meter circuit. A different solution, now under development, entails giving up on the notion of pulse 'counting', *per se*, and replacing the standard threshold discriminator with a new circuit combining the functions of a threshold discriminator, a *pedestal generator*, and a *linear gate* ². The sketch below compares the input-output characteristic of a standard discriminator versus our new circuit.

The output of a standard discriminator circuit is zero for input pulses less than the threshold, and steps to a fixed, pre-determined value for input pulses which exceed the threshold. In the new circuit, the output is again zero for input pulses which are less than the threshold; when the input pulse exceeds the threshold, the output steps, *then linearly follows the amplitude of the input*.


The analog time-averaged (analog rate-meter) output signal from this circuit is proportional to the time-average of energy absorbed (i.e., dose-rate) in the detector probe. The new circuit retains the noise-reducing and drift-reducing advantages of a standard threshold discriminator at low count rates, but with the added advantage that integrated energy/amplitude information contained in

¹G.H. White "The Generation of Random-Time Pulses at an Accurately Known Mean Rate and Having a Nearly Perfect Poisson Distribution" J. Sci Instrum. 1964, Vol 41

² W.R. Leo; Chapter 14.6 in Techniques for Nuclear and Particle Physics Experiments: A How-To Approach, Springer Verlag, ISBN 0-387-57280. New York, Berlin, Heidelberg, 1994

signal pulses which overlap and 'pile up' is preserved over a substantially greater range of input excitations. Our useful range now extends well beyond the point where a standard discriminator's output has 'flat-lined'.

The plots below compare three different detector outputs versus input count rate excitation. The vertical scales are normalized so that all the curves are tangent at low input count rates. In our present system, 'busy time' for a single event is governed by the shaping-amplifier's pulse-width, which is on the order of ~25 micro-seconds – in our case a necessary but reasonable compromise between low dead-time and low noise floor. A wider system bandwidth (shorter shaping time-constant) would allow a narrower pulse which, in turn, would yield a higher maximum count rate, but that would come at the cost of a higher noise floor, requiring a correspondingly higher threshold setting, potentially compromising performance for lower-energy photon-emitters.




As shown below, the *digital* output count-rate peaks at ~17 kHz for 50 kHz input, then gradually declines due to a 'paralyzing dead-time' component³ and finally plateaus at ~13 kHz. However, the analog-rate-meter – or analog average – of that same time-over-threshold discriminator signal has a significantly greater dynamic range, since the discriminator's output pulses vary in duration, staying 'high' when responding to multiple, overlapping input pulses as long as they are of sufficient amplitude to exceed the pre-set threshold. Of course the time-over-threshold analog-rate-meter's output eventually saturates as well, but with a gradual and asymptotic, 'non-paralyzing' characteristic.

New Circuit Our new discriminator circuit significantly extends the useable range of the detector. With this circuit, saturation effects begin to set in at ~150 kHz input count-rate, but the analog output is monotonic – still increasing – up to the present limit of our test apparatus.

The simplest, most common means to achieve detector system DC base-line stability – absolutely vital at low count-rates – is to employ capacitive AC coupling with base-line restoration at the input to the discriminator. That, however, combined with the shaping amplifier's constrained bandwidth, leads to a loss of 'DC-average' information, ultimately causing the apparent signal drop-off at high count rates.

We are currently revisiting many of our prior circuit design assumptions. At the time of this submission, we are seeing preliminary, albeit intriguing and very encouraging test-bench results suggesting there is reason to expect significant improvement over the results posted here.



³ Knoll, Glenn F; Chapter 3, sec. VII in Radiation Detection and Measurement; John Wiley and Sons New York, 1979.